

FREESTYLE SWIM FOOT (LP2-W2)

It's push button easy!

- Users can now walk proudly and comfortably to and from the swim area
- Adjust to 70° plantar flexion for swimming with the simple push of a button
- Utilize a flipper to gain swimming speed and decrease energy expenditure
- Balance the body and enhance swimming efficiency
- Rated to 255 lbs., 12-month warranty, 60 day trial period



FREESTYLE SWIM FOOT
(LP2-W2)

FREESTYLE SWIM FOOT (LP2-W2)

Size: 22-30 cm


User Weight Rating: 255 lbs

Stiffness Category: 1-6

Connector: Male pyramid

Heel Height: 3/8"

Warranty: 12 months



Clearance & Weight:

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	3.99"	485g	127g
26-28 cm	4.07"	542g	159g
29-30 cm	4.55"	635g	219g

Note: Spectra Sock and connector are included with purchase.

REQUIRED ACCESSORIES

Cosmetic Foot Shell:

Mid-profile (FTC-2M)

OPTIONAL ACCESSORIES

- **Shell Removal Tool- Part #ACC-00-103**
- **Additional Spectra Sock- Part #ACC-00-101xx (12", 14", 18")**
- **Additional Stiffening Bumpers-**
 - ACC-00-12201 40 Durometer, Sizes 22-25cm
 - ACC-00-12202 40 Durometer, Sizes 26-28cm
 - ACC-00-12203 40 Durometer, Sizes 29-30cm