


LP SYMES PROSTHETIC FOOT (LP2)

Designed for the demanding needs of a Symes amputee

- Carbon fiber construction delivers energy return and durability
- Remarkably lightweight—averaging 295g
- Direct Lamination Cup allows for less than 2.1" of clearance
- Rated up to 365 lbs



LP SYMES (LP2)

LP SYMES PROSTHETIC FOOT (LP2)

Size*: 22-30 cm

*Special Order Size 31 cm (7-10 day lead time)


User Weight Rating: 365 lbs

Stiffness Category: 1-9

Connector: Direct lamination cup

Heel Height: 3/8"

Warranty: 36 months



Clearance & Weight:

| Size | Clearance | Avg. Weights | |
|----------|-----------|--------------|-------|
| | | Foot | Shell |
| 22-25 cm | 1.7 in | 227g | 127g |
| 26-28 cm | 1.9 in | 295g | 159g |
| 29-30 cm | 2.1 in | 363g | 219g |

Note: A Spectra sock and direct lamination cup are included with purchase.

REQUIRED ACCESSORIES

Cosmetic Foot Shell:

Mid-profile (Default) FTC-2M

High-profile (Optional) FTC-2H

Choice of caucasian or ebony skin tone

OPTIONAL ACCESSORIES

- **Foot Shell Removal Tool- Part #ACC-00-103**
- **Additional Spectra Sock- Part #ACC-00-101xx (12", 14", 18")**
- **Replacement Direct Lamination Cup- Part #ACC-00-162**
- **Additional Stiffening Bumpers-**
 - KIT-00-111 40 Durometer, Sizes 22-25cm
 - KIT-00-113 40 Durometer, Sizes 26-28cm
 - KIT-00-115 40 Durometer, Sizes 29-31cm