

## RENEGADE® PROSTHETIC FOOT (RS1)

*The industry's most innovative shock absorbing prosthetic foot*

- Patented Z-Shock Technology™ achieves unprecedented efficiency, allowing active tibial progression
- Less than one-third the weight of alternative shock foot products
- The all-carbon, z-shaped heel is the most durable shock design available
- Carbon fiber lay-ups are engineered to deliver a smooth gait at any speed, resulting in a comfortable walking foot that can stand up to high-impact forces
- Incorporates full-length sole plate, creating natural inversion-eversion
- Rated up to 365 lbs



RENEGADE® (RS1)

# RENEGADE® (RS1)

**Size\*: 22-30 cm**

\*Special Order Size 31 cm (7-10 day lead time)


**User Weight Rating:** 365 lbs

**Stiffness Category:** 1-9

**Heel Height:** 3/8"

**Connector:** Integral male pyramid

**Warranty:** 36 months



## Clearance & Weight:

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	5.8 in	358g	127g
26-28 cm	6.6 in	437g	159g
29-31 cm	7.3 in	549g	219g

Note: A Spectra sock and pyramid connector are included with the Renegade purchase.

## REQUIRED ACCESSORIES

### Cosmetic Foot Shell:

- Mid-profile (Default) FTC-2M
- High-profile (Optional) FTC-2H
- Choice of caucasian or ebony skin tone

## OPTIONAL ACCESSORIES

### • Pre-Cut Shin Foam (Two part clamshell design)

AFP-00-RS101 For RS1, Sizes 22-25

AFP-00-RS102 For RS1, Sizes 26-28

AFP-00-RS103 For RS1, Sizes 29-31

### • Replacement Renegade Sole Plate- Part #RS1-SP

### • Foot Shell Removal Tool- Part #ACC-00-103

### • Additional Spectra Sock- Part #ACC-00-101xx (12", 14", 18")

### • Additional Stiffening Bumpers-

KIT-00-116 40 Durometer, Sizes 22-25cm

KIT-00-117 40 Durometer, Sizes 26-28cm

KIT-00-118 40 Durometer, Sizes 29-31cm