

RENEGADE® MX (RS1-M1 OR RS6-M1)

Z-Shock Technology™ with enhanced range of motion = superior stability

- Enhanced range of motion easily overcomes changes in surface angles and terrain to provide stability
- Rotation is absorbed at the foot/ankle system rather than the socket, promoting limb health
- 8° of transverse rotation; 8° degrees of inversion/eversion, 6° of plantar flexion, 7.5° of dorsi flexion
- Optional standard height or low profile configuration
- Rated to 285 lbs



Renegade LP-MX

RENEGADE® MX
(RS1-M1 OR RS6-M1)

RENEGADE® MX (RS1-M1 OR RS6-M1)

Part #: Renegade MX: RS1-M1
Renegade LP-MX: RS6-M1

Size*: 22-30 cm

*Special Order Size 31 cm (7-10 day lead time)


User Weight Rating: 285 lbs

Stiffness Category: 1-7

Heel Height: 3/8"

Connector: Integral male pyramid

Warranty: Carbon - 36 months
Rotator - 24 months



Clearance & Weight (MX):

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	7.2 in	672g	127g
26-28 cm	7.25 in	742g	159g
29-31 cm	7.4 in	794g	219g

Clearance & Weight (LP-MX):

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	5.6 in	553g	127g
26-28 cm	5.65 in	646g	159g
29-31 cm	5.7 in	717g	219g

Note: A Spectra sock and pyramid connector are included with the Renegade purchase.

REQUIRED ACCESSORIES

Cosmetic Foot Shell:

Mid-profile (Default for Renegade MX) FTC-2M

High-profile (Default fo Renegade LP-MX) FTC-2H

Choice of caucasian or ebony skin tone

OPTIONAL ACCESSORIES

- Replacement Renegade Sole Plate- Part #RS1-SP
- Foot Shell Removal Tool- Part #ACC-00-103
- Additional Spectra Sock- Part #ACC-00-101xx (12", 14", 18")
- Additional Stiffening Bumpers-
 - KIT-00-116 40 Durometer, Sizes 22-25cm
 - KIT-00-117 40 Durometer, Sizes 26-28cm
 - KIT-00-118 40 Durometer, Sizes 29-31cm