

SILHOUETTE® LP PROSTHETIC FOOT (RS5)

High science meets high style with simplified fitting requirements

- A male pyramid eliminates the need to measure and cut a pylon, simplifying the fitting process
- Slim profile and advanced materials technology facilitate easy and refined cosmesis
- Strong, lightweight carbon-fiber construction allows user to walk farther and faster, using less energy
- Multi-axial rotation delivers an extended range of motion and exceptional stability
- Advanced dual heel system provides a dynamic response for normal ambulation and an override spring for security during higher activity
- Ideal for low to moderate impact K3 amputees
- Rated up to 365 lbs



SILHOUETTE® LP (RS5)

SILHOUETTE® LP PROSTHETIC FOOT (RS5)

Size*: 22-30 cm

*Special Order Size 31 cm (7-10 day lead time)


User Weight Rating: 365 lbs

Connector: 34 mm male pyramid clamp with insert

Stiffness Category: 1-9

Heel Height: 3/8"

Warranty: Carbon - 36 months



Clearance & Weight:

Size	Minimum Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	6.8 in	357g	127g
26-28 cm	7.3 in	387g	159g
29-31 cm	7.7 in	417g	219g

Note: Spectra sock included with purchase.

REQUIRED ACCESSORIES

Cosmetic Foot Shell:

Mid-profile (Default) FTC-2M

High-profile (Optional) FTC-2H

Choice of caucasian or ebony skin tone

OPTIONAL ACCESSORIES

- Replacement Male Pyramid Clamp- Part #ACC-00-196
- Replacement 34mm Tube Clamp Insert- Part #ACC-00-124
- Knee Pyramid Adapter, Male (for knees receiving a 30mm pylon)- Part #ACC-00-195
- Replacement Sole Plate, Tri Durometer- Part #RS4-SP
- Foot Shell Removal Tool- Part #ACC-00-103
- Additional Spectra Sock- Part #ACC-00-101xx (12", 14", 18")
- Additional Stiffening Bumpers-

KIT-00-12101 40 & 80 Durometer, Double Sided Tape, & IFU for Silhouette Sizes 22-25cm

KIT-00-12202 40 & 80 Durometer, Double Sided Tape, & IFU for Silhouette Sizes 26-28cm

KIT-00-12203 40 & 80 Durometer, Double Sided Tape, & IFU for Silhouette Sizes 29-31cm