


SENATOR™ PROSTHETIC FOOT

SOPHISTICATED
SIMPLICITY.



The Senator Prosthetic Foot is a lightweight, economical, energy-returning product that shares composite engineering features with our advanced Freedom Series™ prosthetic feet. It is designed to meet the needs of low to moderately active K3 amputees seeking a prosthetic foot that will allow them to ambulate at varied speeds and participate in daily activities such as bicycling, golfing, hiking, or even jogging.

AN EXCELLENT VALUE AT AN AFFORDABLE PRICE

The Senator Prosthetic Foot combines sophisticated carbon fiber with an integral connector for a lightweight, durable, dynamic-response foot that is simple to fit and finish. It is sophisticated simplicity at its best, delivering substantial energy-return and smooth rollover.

MAJOR BENEFITS

- Lightweight, energy-return design facilitates smooth, easy stride so the user can walk farther, faster using less energy
- Integral male pyramid attachment simplifies the fitting process
- Simulated ankle motion provides a stable base of support throughout stance phase
- Foot is not pre-bonded to the foot shell, allowing access to the foot module and simplifying the use of heel wedges to improve alignment
- Six stiffness categories and low or moderate impact levels enhance customization to ensure ultimate comfort and stability
- Ideal for low to moderately active amputees weighing up to 300 lbs (136 kg), enhancing the range of amputees who can utilize this technology
- Low maintenance, durable design provides peace of mind to the clinician and the user

SENATOR™ PROSTHETIC FOOT



CATEGORY SELECTION CHART

Weight (lbs)	Impact Level		Weight (kg)
	Low	Medium	
100-120	1	2	45-55
121-145	2	3	56-66
146-175	3	4	67-80
176-210	4	5	81-96
211-250	5	6	97-114
251-300	6		115-136

PRODUCT CLEARANCE & WEIGHT

Size (cm)	Clearance	Weight	
		Foot	Shell
22-25	4.5 in	362g	127g
26-28	5.0 in	390g	159g
29-30	5.4 in	418g	219g

SPECIFICATIONS

Size offering - 22-30cm

Carbon fiber foot warranty - 24 months

Foot shell warranty - 6 months

Reimbursement code - L5976



7 Studebaker, Irvine, CA 92618

Phone 949-672-0032

Toll-free 888-818-6777

Fax 949-672-0084

www.freedom-innovations.com

Also distributed through SPS, contact 800-767-7776

© 2006 Freedom Innovations, Inc. All rights reserved.