



# Renegade<sup>®</sup> A·T & Renegade<sup>®</sup> LP-A·T

## All-Terrain Performance AMPLified

- The most dynamic, versatile and durable Renegade prosthetic foot
- EnduraCore<sup>®</sup> Technology delivers 3 times the fatigue life of standard laminates
- Unique design returns up to 35% more energy than competitive alternatives
- Split-toe design delivers multi-axial function with +/- 15 degrees of inversion/eversion for enhanced stability
- Z-Shock Technology<sup>™</sup> absorbs shear forces providing greater user comfort
- Lightweight construction – only 515 grams
- Rated to 166 kg (365 lbs)
- 36-month warranty



Renegade A·T (R11)

Renegade LP-A·T (R16)



Variable Stiffness



Split Toe



Multi-Axial



ACTIVE

**Product ID:** Renegade A-T (R11); Renegade LP-A-T (R16)

**Sizes:** 22-31 cm

**User Weight Rating:** 166 kg (365 lbs)

**Heel Height:** 10 mm (3/8")

**Connector:** Male pyramid

**Warranty:** Foot module 36 months,  
Cosmetic foot shell 6 months


## Build Height & Weight

	Size (cm)	Height (in / mm)	Weight (oz / g)
Renegade A-T	22-25	6 <sup>3</sup> / <sub>8</sub> / 162	22 / 625
	26-28	7 <sup>1</sup> / <sub>4</sub> / 184	27 / 765
	29-31	7 <sup>7</sup> / <sub>8</sub> / 200	33 / 945
Renegade LP-A-T	22-25	5 / 127	22 / 625
	26-28	5 <sup>1</sup> / <sub>2</sub> / 140	27 / 760
	29-31	5 <sup>7</sup> / <sub>8</sub> / 149	33 / 930

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

## Stiffness Category Selection

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9		131-147
326-365	9			148-166



Note: Build height measurements represent a fully assembled unit.

## Accessories

- Pre-Cut Shin Foam: (R11 only; two piece clamshell design)
  - 22-25 cm (AFP-00-RS101)
  - 26-28 cm (AFP-00-RS102)
  - 29-31 cm (AFP-00-RS103)
- Sole Plate: (size and category specific)
  - Universal toe: (RS1-SP)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 16"
- Stiffening Bumpers:
  - 22-25cm: (KIT-00-116)
  - 26-28cm: (KIT-00-117)
  - 29-31cm: (KIT-00-118)

## Foot Module Part Number

-00-0

A

-RU

**Product ID**  
 R11 (Renegade AT)  
 R16 (Renegade LP-A-T)

**Category**  
 1-9

**Size**  
 22-31

## Foot Shell Part Number

FTC-

-

4-R

**Generation**  
 2-3

**Profile**  
 M-mid  
 H-high

**Cap**  
 0=yes  
 1=no

**Size**  
 22-31

**Skin Tone**  
 L-light  
 M-medium  
 D-dark

**Side**  
 L-left  
 R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.