

WalkTek[®]

Redefining K2 stability

- Offers superior stability and support with iMX Technology™
- iMX Technology™ delivers controlled, multi-axial movement and outstanding ground compliance to overcome small community obstacles and changes in surface angles
- A convenient, ankle-height pyramid simplifies the WalkTek fitting process by making adjustments quick and easy
- The WalkTek prosthetic foot features heavy-duty, graphite composite construction and an extremely durable, low-maintenance design, providing you peace of mind
- Rated to 136 kg (300 lbs)
- 12-month warranty



WalkTek (VS3)



Multi-Axial



Sandal Toe



SPECIALTY

Product ID: VS3

Sizes: 21-30 cm


Sandal toe available for sizes 22-28 cm

User Weight Rating: 136 kg (300 lbs)

Heel Height: 10 mm (3/8")

Connector: Male pyramid

Warranty: Foot module 12 months,
Cosmetic foot shell 6 months



Note: Build height measurements represent a fully assembled unit.

Build Height & Weight

Size (cm)	Height (in / mm)	Weight (oz / g)
21-25	3 ³ / ₈ / 85	17 / 490
26-30	3 ³ / ₄ / 97	22 / 630

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Weight (lbs)	Category	Weight (kg)
100-130	1	45-59
131-180	2	60-81
181-240	3	82-109
241-300	4	110-136

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 16"

Foot Module Part Number

VS3-00-0 A -

Category: 1-4 Size: 21-30 Toe Configuration: RU-regular toe, universal; SR-sandal toe, right; SL-sandal toe, left

Foot Shell Part Number

FTC- - 4 -

Generation: 2-3 Profile: M-mid, H-high Cap: 0-yes, 1-no Size: 21-30 Skin Tone: L-light, M-medium, D-dark Toe: R-regular, S-sandal Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.

FEET

