



Catapult[®] Running

Xtreme energy return

- Dual spring technology increases the amount of energy stored and optimizes the timing of its release – providing up to 30% greater energy return
- Energy propulsion allows users to expend less energy in a more stable and balanced design
- Pyramid attachment facilitates alignment changes and the use of a single prosthetic socket for all activities
- Enhanced user customization with interchangeable soft, medium and firm PowerSprings
- Rated up to 166 kg (365 lbs)
- 12-month warranty




Catapult (FX6)



EXTREME SPORTS

- Product ID:** Jogger (FX6)
- User Weight Rating:** 166 kg (365 lbs)
- Connector:** Male pyramid
- Lead Time:** Contact Customer Service
- Warranty:** 12 months
- Return Policy:** 30 days
(foot and standard tread cannot be modified)



Note: Build height measurements represent a fully assembled unit.

Build Height & Weight

Product	Height (in / mm)	Weight (oz / g)
FX6	91/2 / 241	18 / 510

Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Three Power Springs of correlating stiffness categories will be shipped with each Catapult to further enable user optimization.

Weight (lbs)	FX6		Weight (kg)
	Jogger (1)	Runner (2)	
100-115	1	2	44-52
116-130	2	3	53-59
131-150	3	4	60-68
151-170	4	5	69-77
171-195	5	6	78-88
196-220	6	7	89-100
221-255	7	8	101-116
256-285	8	9	117-130
286-325	9	10	131-147
326-365	10		148-166

1) Jogger: Occasional to frequent recreational running

2) Marathoner: Frequent distance running and training for endurance races

Accessories

- Replacement PowerSprings:
 - FX6-PS: (Jogger)
- Tread, Replacement: (FX6-00-P0014-00)

Foot Module Part Number

-00-A00-RU

Product ID
FX6 (Jogger)

Category
01-10

FEET