

Product Reference Chart



	Product	Product Number	Mobility	Articulating Ankle	Vertical Shock	Multi-Axial	Split Keel	Sandal Toe	Heel Height Adjustable	Sizes	Build Height	Weight Rating	Product Weight (including foot shell) (AVG/Range)	Suggested Lcodes
L5968 Ankle	Kinterra®	[ROM]	K ₃	✓			✓	✓		22 – 31cm	4 3/4 – 4 7/8"	275lbs	755 – 999g	L5981 / L5968
	Kintrol™	[VS4]	K ₂	✓				✓		23 – 30cm	4 3/10 – 4 2/5"	300lbs	859g – 1063g	L5972 / L5968
L5987 Shock Feet	Agilix®	[F15]	K ₃		✓	✓	✓	✓		22 – 31cm	5 3/4 – 6 3/4"	365lbs	533 – 772g	L5987 / L5986
	Maverick™ XT	[F11]	K _{3/4}		✓	✓		✓		23 – 30cm	6 2/5 – 7 1/5"	365lbs	646 – 891g	L5987 / L5986
	Maverick™ XT/AT	[F21]	K _{3/4}		✓	✓	✓	✓		23 – 30cm	6 2/5 – 7 1/5"	365lbs	646 – 891g	L5987 / L5986
	Renegade®	[RS1]	K _{3/4}		✓			✓		22 – 31cm	6 3/8 – 7 7/8"	365lbs	625 – 945g	L5987
	Renegade® LP	[RS6]	K _{3/4}		✓			✓		22 – 31cm	5 – 5 7/8"	365lbs	625 – 945g	L5987
	Renegade® AT	[R11]	K _{3/4}		✓	✓	✓			22 – 31cm	6 3/8 – 7 7/8"	365lbs	625 – 945g	L5987 / L5986
	Renegade® LP-AT	[R16]	K _{3/4}		✓	✓	✓			22 – 31cm	5 – 5 7/8"	365lbs	625 – 945g	L5987 / L5986
	Silhouette® VS	[R14]	K ₃		✓	✓				22 – 31cm	7 – 7 7/8"	365lbs	525 – 785g	L5987 / L5986
	Silhouette® LP-VS	[R15]	K ₃		✓	✓				22 – 31cm	6 5/8 – 7 1/2"	365lbs	430 – 685g	L5987 / L5986
L5980 Feet	DynAdapt®	[F10]	K ₃			✓	✓	✓		22 – 31cm	6 1/5 – 7 1/5"	365lbs	535 – 750g	L5980 / L5986
	Silhouette®	[RS4]	K ₃			✓				22 – 31cm	7 – 7 7/8"	365lbs	525 – 785g	L5980 / L5986
	Silhouette® LP	[RS5]	K ₃			✓				22 – 31cm	6 5/8 – 7 1/2"	365lbs	430 – 685g	L5980 / L5986
L5981 Feet	Maverick Comfort AT	[F22]	K ₃			✓	✓	✓		23 – 30cm	4 3/5 – 4 7/10"	365lbs	626 – 841g	L5981 / L5986
	Thrive™	[FS5]	K _{3/4}				✓			22 – 31cm	5 3/4 – 7"	365lbs	630 – 870g	L5981
	Sierra®	[FS1]	K ₃			✓	✓	✓		22 – 31cm	5 3/4 – 7 3/8"	365lbs	515 – 750g	L5981 / L5986
	Highlander® / Max	[FS3 / FS3-H5]	K ₃				✓	✓		22 – 31cm	4 1/4 – 5 1/8"	365lbs / Max 500lbs	470 – 750g	L5981
	Pacifica®	[FS2]	K ₃				✓	✓		24 – 31cm	2 3/8 – 2 1/2"	365lbs	485 – 720g	L5981
	Pacifica® LP	[FS4]	K ₃				✓	✓		22 – 31cm	1 7/8 – 2 3/8"	365lbs	465 – 760g	L5981
	LP Symes™	[LP2]	K _{3/4}							22 – 31cm	1 3/4 – 2 1/8"	365lbs	420 – 670g	L5981
	Runway®	[RS2]	K ₃					✓	✓	22 – 28cm	3 5/8 – 3 7/8"	255lbs	690 – 792g	L5981 / L5990
	Runway HX®	[RS2-K1]	K ₃					✓	✓	22 – 28cm	3 5/8 – 3 7/8"	255lbs	690 – 792g	L5981 / L5990
	Promenade®	[VS2]	K ₃				✓			19 – 31cm	4 5/8 – 5 3/8"	325lbs	490 – 710g	L5981
Defender™	[FS3]	K ₃							16 – 21cm	2 3/8 – 3"	100lbs	175 – 225g	L5981	
L5976 Feet	Senator®	[VS1]	K ₃						22 – 31cm	5 1/8 – 6"	300lbs	520 – 780g	L5976	
L5972 K2 Feet	WalkTek®	[VS3]	K ₂			✓		✓	21 – 30cm	3 3/8 – 3 3/4"	300lbs	490 – 630g	L5972 / L5986	
Specialty Feet	Catapult® Running	Jogger [FX6]	K _{3/4}							N/A	9 1/2"	365lbs	510g	L5999
	Catapult® Running	Sprinter [FX7]	K _{3/4}							N/A	9 1/2"	365lbs	510g	L5999
	Nitro® Tall	[FS6]	K _{3/4}							N/A	9 1/2"	365lbs	275g	L5999
	Nitro® Short	[FS9]	K _{3/4}							N/A	7 1/4"	365lbs	300g	L5999
	Slalom™ Ski	[SKI]	K _{3/4}							custom	5 1/2"	250lbs	1055g	L5999
	Freestyle™ Swim	[LP2-W2]	K _{3/4}							22 – 31cm	4 – 4 5/8"	255lbs	680 – 945g	L5981/L5999
Knees	Plié® 3 Microprocessor Knee	[KS3 / KD3]	K _{3/4}							N/A	9.25" (Pyramid) 8.75" (Threaded)	275 lbs	2.7lbs / 1235g	L5828 / L5845 L5848 / L5856

