



Prosthetic Component Kits

Prosthetic Component Kits

The Freedom Transtibial and Transfemoral Prosthetic Component Kits are convenient and a great value for the Prosthetist. We offer a comprehensive selection of high quality components for multiple build height requirements. Made of titanium with aluminum tubes, our components are the optimal combination of light weight and high performance.



1. Tall Transtibial Kit

Build Height: 4 1/4 – 16 7/8" (11-43cm)

Description	Part Number
Tube & Tube Clamp	
<input type="checkbox"/> 1. Tube w/Titanium Female Pyramid, 30mmx400mm Aluminum	ACC-00-12800-00
<input type="checkbox"/> 2. Tube Clamp, 30mm, Female Pyramid	ACC-00-12500-00
Choose 1 Pyramid	
<input type="checkbox"/> 3. Pyramid, 4 Hole, Male	ACC-00-13000-00
<input type="checkbox"/> 4. Pyramid, 4 Prong, Male	ACC-00-16000-00

Weight Rating: 300lbs

KIT-TT1-12 **000**
 Pyramid
 3. 4 Hole
 4. 4 Prong

2. Medium Transtibial Kit

Build Height: 2 7/8 – 5" (7-12.5cm)

Description	Part Number
Tube & Tube Clamp	
<input type="checkbox"/> 1. Tube w/Female Pyramid, 30mmx100mm Titanium	ACC-00-12700-00
<input type="checkbox"/> 2. Tube Clamp, 30mm, Female Pyramid	ACC-00-12500-00
Choose 1 Pyramid	
<input type="checkbox"/> 3. Pyramid, 4 Hole, Male	ACC-00-13000-00
<input type="checkbox"/> 4. Pyramid, 4 Prong, Male	ACC-00-16000-00

Weight Rating: 300lbs

KIT-TT2-12 **000**
 Pyramid
 3. 4 Hole
 4. 4 Prong



3. Short Transtibial Kit

Build Height: 1 5/8 – 3 1/2" (4-9cm)

Description	Part Number
Choose 1 Double Female Adapter	
<input type="checkbox"/> 1. Adapter, Double Female Pyramid, 32mm	ACC-00-17032-00
<input type="checkbox"/> 2. Adapter, Double Female Pyramid, 45mm	ACC-00-17045-00
<input type="checkbox"/> 3. Adapter, Double Female Pyramid, 60mm	ACC-00-17060-00
<input type="checkbox"/> 4. Adapter, Double Female Pyramid, 75mm	ACC-00-17075-00
Choose 1 Pyramid	
<input type="checkbox"/> 5. Male Pyramid, Rotating, 4 Hole Base	ACC-00-13500-00
<input type="checkbox"/> 6. Pyramid, 4 Hole, Male	ACC-00-13000-00
<input type="checkbox"/> 7. Pyramid, 4 Prong, Male	ACC-00-16000-00

Weight Rating: 300lbs

KIT-TT3- **0000**
 Double Adaptor
 1. 32 mm 3. 60mm
 2. 45 mm 4. 75mm
 Pyramid
 5. Rotatable
 6. 4 Hole
 7. 4 Prong





4. Short Transfemoral Kit

Build Height: 0 – 1" (0-2.5cm)

Description	Part Number
Shin Portion - Tube & Tube Clamp	
<input type="checkbox"/> 1. Tube w/Titanium Female Pyramid, 30mmx400mm Aluminum	ACC-00-12800-00
<input type="checkbox"/> 2. Tube Clamp, 30mm, Female Pyramid	ACC-00-12500-00
Thigh Portion - Choose One Base	
<input type="checkbox"/> 3. Base, 4 Hole, Rotating Pyramid	ACC-00-13700-00
<input type="checkbox"/> 4. Base, 3 Prong, Rotating Pyramid	ACC-00-13800-00
Choose 1 Pyramid Insert	
<input type="checkbox"/> 5. Pyramid, Rotating, Male	ACC-00-13900-00
<input type="checkbox"/> 6. Pyramid, Rotating, Female	ACC-00-14000-00

KIT-TF4-12 **00**

Base Insert

3. 4 Hole 5. Male

4. 3 Prong 6. Female

Weight Rating: 300lbs with 4 Hole,
220lbs with 3 Prong

5. Medium Transfemoral Kit

Build Height: 1¾ – 3½" (4.5-9cm)

Description	Part Number
Shin Portion - Tube & Tube Clamp	
<input type="checkbox"/> 1. Tube w/Titanium Female Pyramid, 30mmx400mm Aluminum	ACC-00-12800-00
<input type="checkbox"/> 2. Tube Clamp, 30mm, Female Pyramid	ACC-00-12500-00
Thigh Portion - Choose One Double Female Adaptor	
<input type="checkbox"/> 3. Adapter, Double Female Pyramid, 32mm	ACC-00-17032-00
<input type="checkbox"/> 4. Adapter, Double Female Pyramid, 45mm	ACC-00-17045-00
<input type="checkbox"/> 5. Adapter, Double Female Pyramid, 60mm	ACC-00-17060-00
<input type="checkbox"/> 6. Adapter, Double Female Pyramid, 75mm	ACC-00-17075-00
Choose 1 Base	
<input type="checkbox"/> 7. Base, 4 Hole, Rotating Pyramid	ACC-00-13700-00
<input type="checkbox"/> 8. Base, 3 Prong, Rotating Pyramid	ACC-00-13800-00
Pyramid Insert	
<input type="checkbox"/> 9. Pyramid, Rotating, Male	ACC-00-13900-00



KIT-TF5-12 **90**

Double Adaptor Base

3. 32 mm 5. 60mm 7. 4 Hole

4. 45 mm 6. 75mm 8. 3 Prong

Weight Rating: 300lbs with 4 Hole,
220lbs with 3 Prong

6. Tall Transfemoral Kit

Build Height: 4¼ – 17" (11-43cm)



Description	Part Number
Shin Portion - Tube & Tube Clamp	
<input type="checkbox"/> 1. Tube w/Titanium Female Pyramid, 30mmx400mm Aluminum	ACC-00-12800-00
<input type="checkbox"/> 2. Tube Clamp, 30mm, Female Pyramid	ACC-00-12500-00
Thigh Portion - Tube	
<input type="checkbox"/> 3. Tube w/Titanium Female Pyramid, 30mmx400mm Aluminum	ACC-00-12800-00
Choose 1 Tube Clamp	
<input type="checkbox"/> 4. Tube Clamp, 30mm, Female Pyramid	ACC-00-12500-00
<input type="checkbox"/> 5. Tube Clamp, 30mm, Male Pyramid	ACC-00-12600-00
Choose 1 Base	
<input type="checkbox"/> 6. Base, 4 Hole, Rotating Pyramid	ACC-00-13700-00
<input type="checkbox"/> 7. Base, 3 Prong, Rotating Pyramid	ACC-00-13800-00
Choose 1 Pyramid Insert	
<input type="checkbox"/> 8. Pyramid, Rotating, Male	ACC-00-13900-00
<input type="checkbox"/> 9. Pyramid, Rotating, Female	ACC-00-14000-00

KIT-TF6-123

Tube Clamp Base Insert

4. Female 6. 4 Hole 8. Male

5. Male 7. 3 Prong 9. Female

Weight Rating: 300lbs with 4 Hole,
220lbs with 3 Prong

Freedom Innovations is solely focused on developing world class lower limb solutions in close collaboration with amputees and prosthetists. These solutions encompass technological innovations, service, training, educational resources and consultative support to help ensure user satisfaction.

Returns, Adjustments and Credits

We are committed to the complete satisfaction of the prosthetist and the amputee. We deliver our component kits with a 60 day "Satisfaction Guarantee". If you are dissatisfied with a Freedom standard production product for any reason, you may return it for a full refund within 60 days of the original invoice date. Component kits can be returned as **complete kits only (no individual component returns allowed)**. All returns require a Return Authorization Number (RA#), which is obtained by calling Freedom Innovations Customer Service. No returns will be accepted without an RA #. A restocking fee of 20% will be applied to orders returned after the 60 day satisfaction guarantee.



Global Headquarters

3 Morgan
Irvine, CA 92618
U.S.A.
toll-free 888.818.6777
tel 949.672.0032
fax 949.672.0084

www.freedom-innovations.com

European Head Office

Jaargetijdenweg 4
7532 SX Enschede
The Netherlands
toll-free 00800-2806 2806
tel +31 (0)53-20 30 300
fax +31 (0)53-20 30 305

www.freedom-innovations.eu