




















# Adapters - General Use

|   | Description                                    | Part Number     | Weight Rating | Material        |
|---|--|-----------------|---------------|-----------------|
|    | Tube Clamp, 30 mm, Female Pyramid              | ACC-00-12500-00 | 300lbs/136kg  | Titanium        |
|    | Tube Clamp, 30 mm, Male Pyramid                | ACC-00-12600-00 | 220lbs/100kg  | Titanium        |
|    | Tube w/Female Pyramid, 30 mm x 100 mm          | ACC-00-12700-00 | 300lbs/136kg  | Titanium        |
|    | Tube w/Titanium Female Pyramid, 30 mm x 400 mm | ACC-00-12800-00 | 300lbs/136kg  | Aluminum        |
|   | Pyramid, 4 Hole, Male                          | ACC-00-13000-00 | 300lbs/136kg  | Titanium        |
|  | Pyramid, 4 Hole, Female                        | ACC-00-13200-00 | 300lbs/136kg  | Titanium        |
|  | Male Pyramid, Rotating, 4 Hole Base            | ACC-00-13500-00 | 300lbs/136kg  | Titanium        |
|  | Female Pyramid, Rotating, 4 Hole Base          | ACC-00-13600-00 | 300lbs/136kg  | Titanium        |
|  | Base, 4 Hole, Rotating Pyramid                 | ACC-00-13700-00 | 300lbs/136kg  | Titanium        |
|  | Base, 3 Prong, Rotating Pyramid                | ACC-00-13800-00 | 300lbs/136kg  | Stainless Steel |
|  | Pyramid, Rotating, Male                        | ACC-00-13900-00 | 300lbs/136kg  | Titanium        |

|   | Description   | Part Number     | Weight Rating | Material           |
|---|---|-----------------|---------------|--------------------|
|    | Pyramid, Rotating, Female   | ACC-00-14000-00 | 300lbs/136kg  | Titanium           |
|    | Tube, 30 mm graphite, 4 mm wall, Aluminum Inserts   | ACC-00-14900-00 | 375 lbs/170kg | Graphite, Aluminum |
|    | Pyramid, 4 Prong, Male  | ACC-00-16000-00 | 300lbs/136kg  | Titanium           |
|    | Adapter, Double Female Pyramid, 32 mm   | ACC-00-17032-00 | 300lbs/136kg  | Titanium           |
|   | Adapter, Double Female Pyramid, 45 mm   | ACC-00-17045-00 | 300lbs/136kg  | Titanium           |
|  | Adapter, Double Female Pyramid, 60 mm   | ACC-00-17060-00 | 300lbs/136kg  | Titanium           |
|  | Adapter, Double Female Pyramid, 75 mm   | ACC-00-17075-00 | 300lbs/136kg  | Titanium           |
|  | Alignment Jig Pyramid, Assembly with 10 mm bolt (for use with category 6-9 foot products) | ACC-00-20100-00 | N/A           | Titanium           |