

# Defender<sup>™</sup> Pediatric

Premium performance similar to the Highlander is now available for pediatric and petite users

- Ensures smooth rollover and exceptional ground compliance in a flexible design
- Rated up to 45 kg (100 lbs), sizes 16-21 cm
- 18-month warranty



Defender Pediatric (FS3)



SPECIALTY

**Product ID:** FS3 (P)

**Sizes:** 16-21 cm


**User Weight Rating:** 45 kg (100 lbs)

**Heel Height:** 10 mm (3/8")

**Connector:** 22 mm pediatric male pyramid  
(adult 30 mm available, call for details)

**Lead Time:** 5-7 days

**Warranty:** Foot module 18 months,  
Cosmetic foot shell 6 months



Note: Build height measurements represent a fully assembled unit.

## Accessories

- Spectra Sock: (ACC-00-101xx) 8" or 10"
- Stiffening Bumpers:
  - 16-21 cm: (KIT-00-12201-01)

## Build Height & Weight

Size (cm)	Height (in / mm)	Weight (oz / g)
16-18	2 <sup>3</sup> / <sub>8</sub> / 61	6 / 175
19-21	3 / 76	8 / 225

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

## Stiffness Category Selection

Weight (lbs)	Category	Weight (kg)
20-30	1P	9-13
31-40	2P	14-18
41-50	3P	19-22
51-60	4P	23-27
61-70	5P	28-31
71-80	6P	32-36
81-100	1 Adult	37-45

## Foot Module Part Number

**FS3-00-0**   -

<b>Category</b>	<b>Size</b>	<b>Toe Configuration</b>
1P-6P 1A	16-21	For 16-18 cm: RR-regular toe, right RL-regular toe, left For 19-21 cm: RU-regular toe, universal

## Foot Shell Part Number

**FTC-**  **M-1**   **4-R**

<b>Generation</b>	<b>Size</b>	<b>Skin Tone</b>	<b>Side</b>
1-(16-18 cm) 2-(19-21 cm)	16-21	L-light M-medium D-dark	L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.

