

Freedom DynAdapt™

Fit, Finish and Roll

- Slim profile, carbon fiber foot system with a slender, anatomic design
- Full length attached heel with no bolts or bushings
- Combination of smooth rollover function, improved stability, and dynamic response makes DynAdapt the best all-around foot for all K3 users
- EnduraCore® Technology delivers up to 3 times the fatigue life of a standard laminate
- Split keel and heel design delivers multi-axial function for enhanced stability
- Lightweight construction – only 635 grams / 22.4 oz (incl. foot shell)
- Rated for 166 kg (365 lbs)
- 36-month warranty



Split Toe



Multi Axial



Sandal Toe

Freedom DynAdapt (F10)

Product ID: F10

Sizes: 22-31 cm


Sandal-toe available for sizes 22-28 cm

User Weight Rating: 166 kg (365 lbs)

Heel Height: 10 mm (3/8")

Connector: Male pyramid

Warranty: Foot module 36 months,
Cosmetic foot shell 6 months



Note: Build height measurements represent a fully assembled unit.

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 16"
- Stiffening Bumpers:
 - 22-25cm: (KIT-00-11373-00)
 - 26-28cm: (KIT-00-11374-00)
 - 29-31cm: (KIT-00-11375-00)

Build Height & Weight

Size (cm)	Height (in / mm)	Weight (oz / g)
22-25	6½ / 158	19 / 535
26-28	6¾ / 170	22 / 635
29-31	7½ / 182	26 / 750

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9	Special Order	131-147
326-365	9	Special Order	Special Order	148-166

Foot Module Part Number

F10-00-0 **A** -

Category: 1-9 Size: 22-31 Toe Configuration: RU-Regular toe, Universal, SR-Sandal toe, Right; SL-Sandal toe, Left

Foot Shell Part Number

FTC- - **4-**

Generation: 2-3 Profile: M-mid, H-high Cap: 0=yes, 1=no Size: 22-31 Skin Tone: L-light, M-medium, D-dark Toe: R-Regular, S-Sandal Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.