






Foot Product Build Height & Weight

	Product	Size (cm)	Build Height ¹		Average Weight ²		
			(in)	(mm)	(oz)	(g)	
 ACTIVE	Renegade® (RS1) & Renegade® A-T (R11)	22-25 26-28 29-31	6 ¾ 7 ¼ 7 ¾	162 184 200	22 27 33	625 765 945	
	Renegade® LP (RS6) & Renegade® LP-A-T (R16)	22-25 26-28 29-31	5 5 ½ 5 ¾	127 140 149	22 27 33	625 760 930	
	Thrive® (FS5)	22-25 26-28 29-31	5 ¾ 6 ½ 7	145 157 175	22 27 31	630 760 870	
 EVERYDAY	Freedom Agilix™ (F15)	22-25 26-28 29-31	5 ¾ 6 ½ 6 ¾	143 161 172	19 24 27	553 609 772	
	Freedom DynAdapt™ (F10)	22-25 26-28 29-31	6 ½ 6 ¾ 7 ½	158 170 182	19 22 26	535 635 750	
	Highlander® (FS3)	22-25 26-28 29-31	4 ¼ 4 ¾ 5 ½	108 121 130	1/ 21 26	4/0 590 750	
	Kinterra™ (ROM)	22-25 26-28 29-31	4 ¾ 4 ¾ 4 ¾	120 121 124	28 32 35	795 897 992	
	Pacifica® (FS2)	24-25 26-28 29-31	2 ¾ 2 ¾ 2 ½	59 61 63	17 22 25	485 610 720	
	Pacifica® LP (FS4)	22-25 26-28 29-31	1 ¾ 2 2 ¾	47 50 60	16 20 27	465 570 760	
	Promenade® (VS2)	19-25 26-28 29-31	4 ½ 4 ¾ 5 ¾	116 124 135	17 21 25	490 595 710	
	Senator® (VS1)	22-25 26-28 29-31	5 ½ 5 ½ 6	130 140 152	18 23 28	520 655 780	
	Sierra® (FS1)	22-25 26-28 29-31	5 ¾ 6 ¾ 7 ¾	146 172 187	18 23 26	515 640 750	
	 STYLE	Runway®/ HX (RS2 & RS2-K1)	22-25 26-28	3 ¾ 3 ¾	92 98	22 26	620 735
		Silhouette® (RS4 & RS8) & Silhouette® VS (R14 & R18)	22-25 26-28 29-31	7 7 ½ 7 ¾	178 191 201	19 23 28	525 645 785
		Silhouette® LP (RS5 & RS9) & Silhouette® LP-VS (R15 & R19)	22-25 26-28 29-31	6 ½ 7 ½ 7 ½	168 181 191	15 19 24	430 550 685
Defender™ Pediatric (FS3)		16-18 19-21	2 ¾ 3	61 76	6 8	175 225	
 SPECIALTY	LP Symes™ (LP2) ³	22-25 26-28 29-31	1 ¾ 1 ¾ 2 ½	43 48 53	15 19 24	420 545 670	
	WalkTek® (VS3)	21-25 26-30	3 ¾ 3 ¾	85 97	17 22	490 630	
	 EXTREME SPORTS	Catapult® Running (FX6/FX7)		9 ½	241	18	510
Freestyle™ Swim (LP2-W2)		22-25 26-28 29-31	4 4 ¼ 4 ¾	102 107 117	24 28 33	680 790 945	
Nitro™ Running (FS6/FS9)		tall (FS6) short (FS9)	9 ½ 7 ¼	241 184	11 10	300 275	
Slalom™ Ski (SKI)		custom	5 ½	140	37	1055	
Foot Shell w/sock (average weights)		22-25 26-28 29-31			7 9 11	195 250 310	

¹ Measurements are representative of fully assembled units including foot shell and protective sock(s) if applicable. Build height is measured from the bottom of the heel to the middle of the pyramid skirt for male adapters, or to the top of female adapters. Build height will vary based on category selected.

² Actual weight varies based on category selected. The reported measurements are average weights.

³ The LP Symes foot is measured to the top of the posterior surface of the lamination cup.