

FREEDOM
MAVERICK™
COMFORT AT

Stride Comfortably!

Danny M. (BK)
Tustin, CA

Introducing the Ultimate Fiberglass Walking Experience



Low Foot Profile

Maverick Comfort AT Fiberglass Foot combines a lower profile foot at an average build height of 4.6" with a split keel so those with clearance challenges can stride with ease wherever they want.



High Energy Return

Maverick Comfort AT offers NO compromise on energy storage and return to provide users with all-day comfort even in a lower profile foot.



Ample Ground Compliance

True to Maverick AT characteristics, Comfort AT offers 32° coronal motion for enhanced ground compliance, so users can walk comfortably everywhere all day long.



Target Users

K3

- K3 level users who prefer an alternative to Carbon Fiber for extra durability and flexibility in a walking foot.
- Any amputees who require the offered build height
- Unilateral and bilateral amputees

Technical Specifications

Sizes: 23 – 30 cm

User Weight Rating: 365 lb (166 kg)

Stiffness Categories: 1 – 9

Build Height: 4.6" / 117 mm (size 26 cm, category 4)

Average Product Weight: 1.53 lb / 693 g

(size 26 cm, category 4, including foot shell)

Toe Option: Regular and Sandal

Split Keel: Yes

Foot Shell Skin Tones: Light, Medium, Dark

Warranty: Foot Module 36 months –

Cosmetic foot shell 6 months

Accessories

- **Foot Shell Removal Tool:** (ACC-00-103)
- **Spectra Socks:** (ACC-00-101xx) 10", 12", 14", 16"
- **Stiffening Bumpers:**
 - 23 – 25 cm: (KIT-11-11433-00)
 - 26 – 28 cm: (KIT-00-11434-00)
 - 29 – 30 cm: (KIT-00-11435-00)

www.freedom-innovations.com
www.maverick-mayhem.com

Global Headquarters

3 Morgan
Irvine, CA 92618
toll free: 888-818-6777
phone: 949-672-0032
fax: 949-672-0084
www.freedom-innovations.com



#freedominnovations

European Head Office

Jaargetijdenweg 4
7532 SX Enschede
The Netherlands
toll free: 00800 - 2806 2806
phone: 31 (0) 53 20 30 300
fax: 31 (0) 53 20 30 305
www.freedom-innovations.eu

USA Manufacturing & Returns

425 East 400 North
Gunnison, UT 84634



1 Ordering Guide

Verify User's Specifications

- Patient Weight
- Foot Size (23 -30 cm)
- Right or Left
- Regular Toe or Sandal Toe*

*Available Sandal Toe Sizes: 23 -28 cm

2 Stiffness Category Selection Chart

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9		131-147
326-365	9			148-166

3 Complete Foot Module Part Number

