

FREEDOM  
**MAVERICK™**  
COMFORT AT

Stride Comfortably!

Danny M. (BK)  
Tustin, CA

*Introducing the Ultimate Fiberglass Walking Experience*



#### Low Foot Profile

Maverick Comfort AT Fiberglass Foot combines a lower profile foot at an average build height of 4.6" with a split keel so those with clearance challenges can stride with ease wherever they want.



#### High Energy Return

Maverick Comfort AT offers NO compromise on energy storage and return to provide users with all-day comfort even in a lower profile foot.



#### Ample Ground Compliance

True to Maverick AT characteristics, Comfort AT offers 32° coronal motion for enhanced ground compliance, so users can walk comfortably everywhere all day long.



## Target Users

**K3**

- K3 level users who prefer an alternative to Carbon Fiber for extra durability and flexibility in a walking foot.
- Any amputees who require the offered build height
- Unilateral and bilateral amputees

## Technical Specifications

**Sizes:** 23 – 30 cm

**User Weight Rating:** 365 lb (166 kg)

**Stiffness Categories:** 1 – 9

**Build Height:** 4.6" / 117 mm (size 26 cm, category 4)

**Average Product Weight:** 1.53 lb / 693 g

(size 26 cm, category 4, including foot shell)

**Toe Option:** Regular and Sandal

**Split Keel:** Yes

**Foot Shell Skin Tones:** Light, Medium, Dark

**Warranty:** Foot Module 36 months –

Cosmetic foot shell 6 months

## Accessories

- **Foot Shell Removal Tool:** (ACC-00-103)
- **Spectra Socks:** (ACC-00-101xx) 10", 12", 14", 16"
- **Stiffening Bumpers:**
  - 23 – 25 cm: (KIT-11-11433-00)
  - 26 – 28 cm: (KIT-00-11434-00)
  - 29 – 30 cm: (KIT-00-11435-00)

[www.freedom-innovations.com](http://www.freedom-innovations.com)  
[www.maverick-mayhem.com](http://www.maverick-mayhem.com)

### Global Headquarters

3 Morgan  
Irvine, CA 92618  
toll free: 888-818-6777  
phone: 949-672-0032  
fax: 949-672-0084  
[www.freedom-innovations.com](http://www.freedom-innovations.com)



#freedominnovations

### European Head Office

Jaargetijdenweg 4  
7532 SX Enschede  
The Netherlands  
toll free: 00800 - 2806 2806  
phone: 31 (0) 53 20 30 300  
fax: 31 (0) 53 20 30 305  
[www.freedom-innovations.eu](http://www.freedom-innovations.eu)

### USA Manufacturing & Returns

425 East 400 North  
Gunnison, UT 84634



## 1 Ordering Guide

### Verify User's Specifications

- Patient Weight
- Foot Size (23 -30 cm)
- Right or Left
- Regular Toe or Sandal Toe\*

\*Available Sandal Toe Sizes: 23 -28 cm

## 2 Stiffness Category Selection Chart

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9		131-147
326-365	9			148-166

## 3 Complete Foot Module Part Number

