



1 2 3 4



SPECIALTY

Defender™ Pediatric

Premium performance similar to the Highlander is now available for pediatric and petite users

- Ensures smooth rollover and exceptional ground compliance in a flexible design
- Rated up to 45 kg, sizes 16-21 cm
- 18-month warranty




Defender Pediatric (FS3)

Defender™ Pediatric



SPECIALTY



See Build Height Chart

10 mm

Note: Build height measurements represent a fully assembled unit.

Product ID: FS3

Sizes: 16-21 cm

User Weight Rating: 45 kg

Heel Height: 10 mm

Connector: 22 mm pediatric male pyramid (adult 30 mm available, call for details)

Lead Time: 5-7 days

Warranty: Foot module 18 months, Cosmetic foot shell 6 months

Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
16-18	61	175
19-21	76	225

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Category	Weight (kg)
1P	9-13
2P	14-18
3P	19-22
4P	23-27
5P	28-31
6P	32-36
1 Adult	37-45

Accessories

- Spectra Sock: (ACC-00-101xx)
20 cm (08) or 25 cm (10)
- Stiffening Bumpers:
- 16-21 cm (KIT-00-12201-01)

Note: Spectra socks and heel bumpers are included in the delivery.

Complete Prosthetic Foot (incl. foot shell) Part Number

FS3-00-1 **-M**

Category: 1P-6P, 1A; Size: 16-21; Side: L-left, R-right; Generation: 1-(16-18cm), 2-(19-21cm); Skin Tone: L-light, M-medium, D-dark

Foot Module Part Number

FS3-00-0 **-**

Category: 1P-6P, 1A; Size: 16-21; Toe Configuration: For 16-18 cm: RR-regular toe, right; RL-regular toe, left; For 19-21 cm: RU-regular toe, universal

Foot Shell Part Number

FTC- **M-1** **4-R**

Generation: 1-(16-18 cm), 2-(19-21 cm); Size: 16-21; Skin Tone: L-light, M-medium, D-dark; Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.