

EVERYDAY

Promenade[®]

Designed specifically for low impact amputees

- Near effortless roll-over optimized for low impact amputees, so they can do more and feel better at the end of the day
- Partial split toe offers controlled inversion/eversion, providing users enhanced confidence for a variety of walking conditions
- Full length sole plate facilitates a smooth transition from heel strike all the way through toe off, providing a more fluid gait and increased comfort throughout the day
- Male pyramid and only 5 stiffness categories simplify the fitting process, saving you time and money
- The 147 kg weight rating makes the product appropriate for a broader range of users compared to competitive alternatives
- 36-month warranty




Promenade (VS2)



Split Toe

Promenade®



Note: Build height measurements represent a fully assembled unit.

Product ID: VS2

Sizes: 19-31 cm

User Weight Rating: 147 kg

Heel Height: 10 mm

Connector: Male pyramid

Warranty: Foot module 36 months,
Cosmetic foot shell 6 months

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
 - 22-25 cm (KIT-00-11393)
 - 26-28 cm (KIT-00-11394)
 - 29-31 cm (KIT-00-11395)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
19-25	116	490
26-28	124	595
29-31	135	710

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Impact Level	Weight (kg)
Low	
1	45-57
2	58-72
3	73-92
4	93-117
5	118-147

Complete Prosthetic Foot (incl. foot shell) Part Number

VS2-00-0 A -M -3

Category 1-5 Size 19-31 Side L-left R-right Skin Tone L-light M-medium D-dark

Foot Module Part Number

VS2-00-0 A -RU

Category 1-5 Size 19-31

Foot Shell Part Number

FTC- - 4-R

Generation 2-3 Profile M-mid H-high Cap 0-yes 1-no Size 19-31 Skin Tone L-light M-medium D-dark Side L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.