



1

2

3

4



EVERYDAY

# WalkTek®

## Redefining K1 and K2 stability

- Offers superior stability and support with iMX Technology™
- iMX Technology™ delivers controlled, multi-axial movement and outstanding ground compliance to overcome small community obstacles and changes in surface angles
- A convenient, ankle-height pyramid simplifies the WalkTek fitting process by making adjustments quick and easy
- The WalkTek prosthetic foot features heavy-duty, graphite composite construction and an extremely durable, low-maintenance design, providing you peace of mind
- Rated to 136 kg
- 12-month warranty




WalkTek (VS3)



Multi-Axial



Sandal Toe



Note: Build height measurements represent a fully assembled unit.

**Product ID:** VS3

**Sizes:** 21-30 cm

Sandal toe available for sizes 22-28 cm\*

**User Weight Rating:** 136 kg

**Heel Height:** 10 mm

**Connector:** Male pyramid

**Warranty:** Foot module 12 months,

Cosmetic foot shell 6 months

\*When ordering the Sandal Toe option, please contact Customer Service for lead time.

**Accessories**

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)  
30 cm, 35 cm, 40 cm

Note: Spectra socks and heel bumpers are included in the delivery.

**Build Height & Weight**

Size (cm)	Height (mm)	Weight (g)
21-25	85	490
26-30	97	630

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

**Stiffness Category Selection**

Impact Level	Weight (kg)
1	45-59
2	60-81
3	82-109
4	110-136

**Complete Prosthetic Foot (incl. foot shell) Part Number**

**VS3-00-A--**

Toe: 0-regular, 1-sandal  
 Category: 1-4  
 Size: 21-30  
 Version: M-mid profile, S-sandal toe  
 Side: L-left, R-right  
 Generation: 3-regular toe, 2-sandal toe  
 Skin Tone: L-light, M-medium, D-dark

**Foot Module Part Number**

**VS3-00-0A-**

Category: 1-4  
 Size: 21-30  
 Toe Configuration: RU-regular toe, universal, SR-sandal toe, right, SL-sandal toe, left

**Foot Shell Part Number**

**FTC---4-**

Generation: 1-3  
 Profile: M-mid, H-high, S-sandal toe  
 Cap: 0-yes, 1-no  
 Size: 21-30  
 Skin Tone: L-light, M-medium, D-dark  
 Toe: R-regular, S-sandal  
 Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.