

# Freestyle<sup>™</sup> Swim

## It's push button easy!

- Users can now walk proudly and comfortably to and from the swim area
- Adjust to 70° plantarflexion for swimming with the simple push of a button
- Utilize a flipper to gain swimming speed and decrease energy expenditure
- Balance the body and enhance swimming efficiency
- Rated to 116 kg (255 lbs)
- Warranty: 12-months ankle, 36 months foot module



Freestyle Swim (LP2-W2)



EXTREME SPORTS

**Product ID:** LP2-W2

**Sizes:** 22-31 cm

**User Weight Rating:** 116 kg (255 lbs)

**Heel Height:** 10 mm (3/8")

**Connector:** Male pyramid

**Warranty:** Foot module 36 months, Ankle 12 months,  
Cosmetic foot shell 6 months


## Build Height & Weight

Size (cm)	Height (in / mm)	Weight (oz / g)
22-25	4 / 102	24 / 680
26-28	4¼ / 107	28 / 790
29-31	4⅝ / 117	33 / 945

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

## Stiffness Category Selection

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6		89-100
221-255	6			101-116



Note: Build height measurements represent a fully assembled unit.

## Accessories

- Protective Cover for Ankle: (KIT-00-200)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 16"
- Stiffening Bumpers:
  - 22-25 cm (KIT-00-11100-00)
  - 26-28 cm (KIT-00-11300-00)
  - 29-31 cm (KIT-00-11500-00)

## Foot Module Part Number

LP2-W2-0  A  -RU

Category: 1-6      Size: 22-31

## Foot Shell Part Number

FTC-3M-1   4-R

Size: 22-31      Skin Tone: L-light, M-medium, D-dark      Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.

FEET