

SPORT

Freestyle™ Swim

It's push button easy!

- Users can now walk proudly and comfortably to and from the swim area
- Adjust to 70° plantar flexion for swimming with the simple push of a button
- Utilize a flipper to gain swimming speed and decrease energy expenditure
- Balance the body and enhance swimming efficiency
- Rated to 116 kg
- 12-month warranty



Freestyle Swim (LP2-W2)




EnduraCore®
Technology

Freestyle™ Swim



SPORT



Note: Build height measurements represent a fully assembled unit.

Product ID: LP2-W2

Sizes: 22-31 cm

User Weight Rating: 116 kg

Heel Height: 10 mm

Connector: Male pyramid

Warranty: Foot module 36 months,
Ankle 12 months, Cosmetic foot shell 6 months
Freestyle Swim can be delivered with a special
shaped cover, please contact Customer Service.

Accessories

- Protective Cover for Ankle: (KIT-00-200)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
 - 22-25cm (KIT-00-11100-00)
 - 26-28cm (KIT-00-11300-00)
 - 29-31cm (KIT-00-11500-00)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
22-25	87	690
26-28	103	860
29-31	111	980

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Impact Level			Weight (kg)
Low	Moderate	High	
1	1	2	44-52
1	2	3	53-59
2	3	4	60-68
3	4	5	69-77
4	5	6	78-88
5	6		89-100
6			101-116

Complete Prosthetic Foot (incl. foot shell) Part Number

LP2-W2-1 A -M -3

Category 1-6 Size 22-31 Side L-left R-right Skin Tone L-light M-medium D-dark

Foot Module Part Number

LP2-W2-0 A -RU

Category 1-6 Size 22-31

Foot Shell Part Number

FTC-3M-1 4-R

Size 22-31 Skin Tone L-light M-medium D-dark Side L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.