



1

2

3

4



SPECIALTY

LP Symes[®]

Designed for the demanding needs of a Symes amputee


- Engineered to withstand the high forces of a Symes amputee, while remaining flexible enough to offer superior comfort
- Remarkably lightweight—averaging 295 g
- Direct lamination cup keeps build height low
- Rated up to 166 kg
- 36-month warranty



LP Symes (LP2)



EnduraCore[®]
Technology



Note: Build height measurements represent a fully assembled unit.

Product ID: LP2

Sizes: 22-31 cm

User Weight Rating: 166 kg

Heel Height: 10 mm

Connector: Direct lamination cup

Warranty: Foot module 36 months,
Cosmetic foot shell 6 months

Accessories

- Direct Lamination Cup: (ACC-00-162)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
 - 22-25 cm (KIT-00-111)
 - 26-28 cm (KIT-00-113)
 - 29-31 cm (KIT-00-115)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
22-25	62	440
26-28	68	590
29-31	76	690

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Impact Level			Weight (kg)
Low	Moderate	High	
1	1	2	44-52
1	2	3	53-59
2	3	4	60-68
3	4	5	69-77
4	5	6	78-88
5	6	7	89-100
6	7	8	101-116
7	8	9	117-130
8	9		131-147
9			148-166

Complete Prosthetic Foot (incl. foot shell)
Part Number

LP2-00-1 A -M -3

Category 1-9 Size 22-31 Side L-left R-right Skin Tone L-light M-medium D-dark

Foot Module Part Number

LP2-00-0 A -RU

Category 1-9 Size 22-31

Foot Shell Part Number

FTC-3M-1 4-R

Size 22-31 Skin Tone L-light M-medium D-dark Side L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.