

Freestyle Swim™

It's push button easy!

- Users can now walk proudly and comfortably to and from the swim area
- Adjust to 70° plantarflexion for swimming with the simple push of a button
- Utilize a flipper to gain swimming speed and decrease energy expenditure
- Balance the body and enhance swimming efficiency
- Rated to 116 kg (255 lbs)
- Warranty: 12-months ankle, 36 months foot module



Freestyle Swim (LP2-W2)



EXTREME SPORTS

Product ID: LP2-W2


Sizes:22-31 cm

User Weight Rating: 116 kg (255 lbs)

Heel Height:10mm (3/8")

Connector: Male pyramid

Warranty: Foot module 36 months, Ankle 12 months, Cosmetic foot shell 6 months



Note: Build height measurements represent a fully assembled unit.

Build Height & Weight

Size (cm)	Height (in/ mm)	Weight (oz/g)
22-25	3.4 / 87	24.3 / 690
26-28	4 / 103	30 / 860
29-31	4.4 / 111	34.6 / 980

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6		89-100
221-255	6			101-116

Accessories

- Protective Cover for Ankle: (KIT-00-200)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 16"
- Stiffening Bumpers:
 - 22-25 cm (KIT-00-11100-00)
 - 26-28 cm (KIT-00-11300-00)
 - 29-31 cm (KIT-00-11500-00)

Foot Module Part Number

LP2-W2-0 A -RU

Category 1-6 Size 22-31

Foot Shell Part Number

FTC-3M-1 4-R

Size 22-31 Skin Tone L-light M-medium D-dark Side L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.

FEET